

[amgmicro.com](http://amgmicro.com)



**EVOS<sup>®</sup>** *fi*

**QUICK START GUIDE**

Digital Inverted Microscope  
for **Fluorescence** and **Transmitted**  
**Light** Applications




ZP-PKGA-0494 REV B

## Overview

The **EVOS fl** microscope has two types of controls: mechanical and software. Mechanical controls include the stage X-Y axis knobs, focusing knobs, objective selection wheel, phase annuli selector, and the light cube selection lever. Software controls are located in the **control bar** at the bottom of the display screen. The **channel bar** at the top of the display screen shows the selected light cube or transmitted light position.

**CAUTION**  
**UV Light!**  
DO NOT Look Directly At Light

Avoid exposure to beam and use protective shields.  
**NEVER look directly at UV light!**



## Basic Operation

1. Turn on the microscope with the **power switch** ① on the right side of the base.
2. Plug a USB flash drive into one of the **USB ports** ② on the right side of the microscope.
3. Place the **sample** ③ on the stage, using a vessel holder if needed.
4. Set magnification with the **objective selection wheel** ④ on the front of the microscope.
5. Pull the **light cube selection lever** ⑤ (left side of base) all the way toward the front of the microscope (the Channel Bar will highlight the "Transmitted" position).



2

6. Turn on illumination with the **LIGHT ON button** located on the left side of the control bar.
7. Focus the sample with **focusing knobs**.
8. **Optional:** To take a picture of the transmitted light image, click the **Capture button** on the control bar.
9. Place the **light shield box** on the stage, over the sample.
10. Move the **light cube selection lever** to the desired fluorescence channel (*the Channel Bar will highlight the selected light cube*).
11. With the **Find & Focus tab** active, turn on fluorescence illumination using the **LIGHT ON button**.
12. Adjust the focus as necessary.
13. Adjust the **Illumination Intensity slider** on the control bar as needed.
14. Click the **Capture button**.
15. Repeat steps 10–14 to acquire each fluorescence channel.
16. Click the **Overlay tab** to show all channels in color overlay mode.
17. Adjust **Brightness and Contrast** for each channel to bring them into desired balance.
18. Click the **Save button** to save the color image (*refer to EVOS fl User Guide*).



## Helpful Tips


- ▶ In Find & Focus or Actual Mode, the **Color option** (for EVOS fl) can be turned off to display a grayscale image. This often shows more details than a color image.



For EVOS flc use the **Color Adjustment button** to fine tune your live image Brightness, Contrast, Saturation, and Hue prior to capture.

- ▶ In Find & Focus Mode, the exposure time is set to 100 ms to assist real-time focusing, moving the stage, etc. The illumination level is approximately 60% of the amount used for image capture, in order to minimize photobleaching and phototoxicity. Clicking **Capture** results in brighter illumination and a longer exposure time during image capture to provide a high-quality image.
- ▶ In Actual Mode, turning on the illumination results in full-powered illumination and actual exposure times for live viewing of the sample. With longer exposure times (*more than 200 ms*) there will be a lag between moving the focus knob and seeing the focus change onscreen.

**Note:** The EVOS fl User Guide is on the USB flash drive. You can also download it from the EVOS fl product page at [www.amgmicro.com](http://www.amgmicro.com).



1. Channel bar (*active channel highlighted*)
2. Mode tabs (*Find & Focus, Actual, Overlay*)
3. LIGHT ON/OFF button
4. Illumination Intensity slider
5. Capture button
6. Color option or Color Adjustment button
7. Settings tab
8. Save button
9. Information bar
10. Selected objective



